

AD-A063 238

ARIZONA UNIV TUCSON DEPT OF CHEMISTRY

F/G 7/4

NOVEL MICROCOMPUTER-CONTROLLED POTENTIOMETRIC ANALYSIS SYSTEM.(U)

JAN 79 C MARTIN, H FREISER


N00014-78-C-0001

UNCLASSIFIED

TR-1

NL

OF
AD
A063238



END

DATE

FILMED

3-79

DDC

~~SECRET~~ 6 LEVEL II

AD A063238

DDC FILE COPY.

OFFICE OF NAVAL RESEARCH

Contract ¹⁵ N00014-78-C-0001

Task No. NR 051-672

TECHNICAL REPORT NO. 1

6 NOVEL MICROCOMPUTER-CONTROLLED
POTENTIOMETRIC ANALYSIS SYSTEM

by

10 Charles/Martin ~~and~~ Henry/Freiser

14 TR-1

Prepared for Publication
in

Analytical Chemistry

9 Technical rept, 1 Apr-31 Dec 78,

DDC
RECEIVED
JAN 16 1979

University of Arizona
Department of Chemistry
Tucson, Arizona 85721

11 5 Jan ~~1979~~ 1979

12 26p

Reproduction in whole or in part is permitted for
any purpose of the United States Government
This document has been approved for public release
and sale; its distribution is unlimited

79 01 12 046

033 860

JOB

REPORT DOCUMENTATION PAGE		READ INSTRUCTIONS BEFORE COMPLETING FORM
1. REPORT NUMBER Number 1	2. GOVT ACCESSION NO.	3. RECIPIENT'S CATALOG NUMBER
4. TITLE (and Subtitle) Novel Microcomputer-Controlled Potentiometric Analysis System		5. TYPE OF REPORT & PERIOD COVERED Technical Report 1 April 78 - 31 Dec 78
7. AUTHOR(s) Charles Martin and Henry Freiser		6. PERFORMING ORG. REPORT NUMBER
9. PERFORMING ORGANIZATION NAME AND ADDRESS University of Arizona Tucson, AZ 85721		8. CONTRACT OR GRANT NUMBER(s) N00014-78-C-0001
11. CONTROLLING OFFICE NAME AND ADDRESS Office of Naval Research		10. PROGRAM ELEMENT, PROJECT, TASK AREA & WORK UNIT NUMBERS
14. MONITORING AGENCY NAME & ADDRESS (if different from Controlling Office)		12. REPORT DATE 5 January 1979
		13. NUMBER OF PAGES 22
		15. SECURITY CLASS. (of this report)
		15a. DECLASSIFICATION/DOWNGRADING SCHEDULE
16. DISTRIBUTION STATEMENT (of this Report) Approved for Public Release; Distribution Unlimited		
17. DISTRIBUTION STATEMENT (of the abstract entered in Block 20, if different from Report) Same		
19. SUPPLEMENTARY NOTES None		
21. KEY WORDS (Continue on reverse side if necessary and identify by block number) Microcomputer Automated Potentiometric Analysis Multiple Ion Selective Electrodes Automated Titration		
23. ABSTRACT (Continue on reverse side if necessary and identify by block number) A novel microcomputer-controlled potentiometric analysis system is described. This system combines the low cost and ease of construction advantages of microcomputer-based instrumentation with the high level software capabilities of a minicomputer. This is made possible through the use of a powerful new computer language called CONVERS. The instrumentation and CONVERS software system are described. Applications of this system to endpoint determination in potentiometric titrations and fully automated generation of calibration curve data for		

ACCESSION for	
NTIS	White Section <input type="checkbox"/>
DOC	Buff Section <input type="checkbox"/>
UNANNOUNCED <input type="checkbox"/>	
JUSTIFICATION _____	
BY _____	
DISTRIBUTION/AVAILABILITY CODES	
Dist. _____	SPECIAL _____
A	

79 01 12 046

INTRODUCTION

A number of examples of computer automation of instrumentation used with ion selective electrodes (ISEs) have appeared in the literature (1-9). The computer-based systems described have ranged from a relatively inflexible single-task potentiometric titration system (3), to an elegant multi-task system allowing a very high degree of interaction between computer and user (5). From the standpoint of either the routine user of ISEs, who is interested in automating analyses currently done manually, or the researcher in the field of ISE development who is interested in efficient and convenient evaluation of either new ion selective sensors or new methodology, nearly all of the systems described, thus far, suffer from one or both of the following disadvantages.

The first disadvantage is expense. Systems that have had sufficient computation ability to accomplish, for example, Gran type (5, 6, 7) or standard addition (4, 8) analysis, have been based on minicomputers which are still very expensive (\$8,000 to \$20,000). For this reason, incorporation of a minicomputer into a dedicated potentiometric analysis system is undesirable. If a laboratory minicomputer is already available, mini-micro combinations (2) are an attractive alternative.

The second disadvantage is complexity. The average ISE user may not possess the electronic and computer expertise to assemble a system capable of meeting analytical demands. This problem is, once again, largely associated with the use of a minicomputer in that due to the total system approach to minicomputer construction, interfacing the computer to external devices is a relatively complex process. This problem is frequently made even more difficult by the need to extend existing input/output (I/O) capabilities through multiplexing.

Microcomputers offer distinct advantages over minicomputers with regard to both cost and complexity. Microcomputers and associated supporting hardware are quite inexpensive, making incorporation into dedicated instrumentation very desirable. Furthermore, because of the modular approach to microcomputer architecture, construction of an analytical system based on a microcomputer is surprisingly easy. Traditionally, microcomputers suffer from lack of software support, however, making high level mathematical manipulations and computer-user interaction very difficult. For this reason, instrumentation based on stand-alone microcomputers described thus far have lacked power and versatility (1, 3).

This paper describes a novel microcomputer-based potentiometric analysis system that incorporates the cost and simplicity advantages of the microcomputer with the high level software capabilities associated with minicomputers. This is made possible through the use of a new computer language, CONVERS (10). Illustrative applications of this system to endpoint determination in potentiometric titrations and fully automated generation of ISE calibration data are presented.

EXPERIMENTAL

Electronic Measurements. An Orion 801A digital pH-millivolt meter was used to measure electrode potentials. When in the pH mode, this meter can measure from pH 0.000 to pH 13.999 in 0.001 pH unit increments. When in the millivolt (mv) mode, it can measure potentials from +1999.9 mv to -1999.9 mv in 0.1 mv increments. Interfacing this device to the computer is accomplished through a convenient TTL-compatible interface tab.

The meter is used in conjunction with an electrode multiplexer similar to that described by Zipper, et. al. (8). This device can switch up to six

different electrodes to the meter either manually via a set of front panel switches or under computer control.

Since Orion does not supply the user with information regarding the analog to digital conversion (ADC) frequency or the noise rejection capability of the 801 meter, these parameters were determined experimentally. The ADC frequency was determined by observing the print gate frequency (interface tab pin 16) on an oscilloscope. The noise rejection capability was evaluated by determining the standard deviation of the digitized signal when a constant (but undoubtedly noisy) analog signal was presented to the meter. An Orion fluoride (94-09) ISE in conjunction with a double junction (Thomas 4092-H10) reference electrode (0.1M KCl outer solution) was used to generate the analog signal. The measuring solution was 0.5M NaF maintained at $25^{\circ}\text{C} \pm 0.1^{\circ}\text{C}$. The electrodes were immersed in this solution one hour prior to measurement to ensure that a constant EMF was being produced. The measuring solution was covered to prevent evaporative losses.

Standard deviations for both single and averaged ADCs were calculated. Simple CONVERS computer algorithms were used to generate the data and perform the necessary calculations. In the single ADC algorithm, the computer holds the meter, samples, frees the meter, waits three seconds and repeats the process. Twenty readings are taken and used to calculate a standard deviation. Ten standard deviations are calculated and averaged during an analysis. In the averaging algorithm, the computer obtains a number of ADCs at the maximum ADC frequency. This is done by continuously checking the print gate and sampling once after each low-high transition (low = busy state). These ADCs are then averaged to give a reading. As before, twenty readings are used to calculate a standard deviation and ten standard deviations are calculated and averaged per analysis.

Buret System. A Mettler DV10 digital buret drive with 10 ml buret (DV101) was used. The accuracy of the buret was determined gravimetrically with water. The buret was found to be accurate to about 0.1% over the entire volume range. The buret tip was submerged in the measuring solution throughout both the titrimetric and calibration experiments. The buret is triggered by a pulse train from the computer. Each pulse causes the buret to deliver 1 μ L of titrant. The maximum pulsing frequency is 166 Hz. An interface board output strobe was used in conjunction with a 74121 monostable integrated circuit chip to provide the required pulse train. The magnetic stirrer is turned on and off by a TTL triggerable relay circuit. The relay switches 120V AC to the stirrer.

Computer and Associated Hardware. A schematic diagram of the computer and associated hardware is shown in Figure 1. All components except the electrode multiplexer are commercially available. The computer is an IMSAI 8080 microcomputer kit (IMSAI Manufacturing, San Leandro, CA) equipped with 24K of static memory (Jade Computer Products, Hawthorne, CA). All interfacing is done through Processor Technology (Emeryville, CA) 3P+S I/O boards which fit directly into the computer's S-100 motherboard. Two 3P+S boards are used. The floppy disk drive is a North Star (Berkeley, CA) MDS-A. The disk drive control board is also S-100 compatible. The X-Y recorder is a Hewlett-Packard 7001A.

Programming. All programming was done in a recently developed computer language, CONVERS, (10). This language is based on an interpretive compiler and combines the conversational and easy debugging characteristics of an interpreter type language (for example, BASIC) with the speed associated with a compiler type language (such as FORTRAN). In addition, I/O control and high level arithmetic routines are provided. The high level I/O routines are a particularly attractive feature of this language in that they obviate the need

for assembly language subroutines for I/O operations.

To demonstrate the capabilities of this system, software to locate endpoints in a potentiometric titration and plot out the titration curve and software to calibrate up to six ISEs and plot out the calibration curves was written. Programs were written with the aid of the CONVERS disk operating system (CDOS) and text editor. The CDOS was used for program storage and retrieval. Overlays were used to conserve memory space. For most applications described, 16K of memory was sufficient.

The potentiometric titration program uses routines similar to that described by Christiansen, et. al. (3) to calculate desired titrant increment volume (DTIV) and locate endpoints. DTIV is calculated on the basis of the magnitude of the first derivative of the titration curve. The endpoint is located by calculating the second derivative and extrapolating between the last positive and first negative values. Two slightly different routines were used, one for acid-base titrations in which the meter is used in the pH mode and one for general potentiometric titrations in which the mv mode is used. The software quizzes the meter to determine which mode is operative and implements the appropriate routine.

A very high degree of computer-user interaction is easily obtained with the CONVERS language. This is evidenced in Figure 2, which shows a typical computer printout from a potentiometric titration (titration of a mixture of acetic and hydrochloric acids with sodium hydroxide). The statements before the equal signs are typed by the computer, the numbers are entered by the user. The concentration entered is in molarity and the sample volume in mL. The points the computer is instructed to take after an endpoint is located allows for detection of multiple endpoints. The titration will cease after the desired number of extra increments are added unless another region of increasing

slope is detected in which case the desired extra increments are added after this additional endpoint has been located.

The drift parameter is entered in mv. The computer will continue sampling the electrodes after an increment has been added until two consecutive readings (separated by 6 seconds) differ by less than this specified amount or until the number of 6 second intervals entered in the next statement has been surpassed. In this case, a maximum drift rate of 0.1 mv per minute and a 1-minute equilibration time were specified. If the electrode drift exceeds these criteria after addition of a particular increment, the computer informs the user and continues with the titration. As indicated in Figure 2, excessive drift was present during this titration after addition of increment number 51.

The maximum and minimum pulses determine the maximum and minimum increments the computer is to add. The maximum increment is added when the algorithm used to calculate increments produces a number greater than the maximum specified and the minimum increment is added when a number smaller than the specified minimum is calculated. For this particular titration, the maximum increment was rather small (70 μ l) so that the poorly defined first endpoint corresponding to the concentration of HCl, could be located (see Figure 3). Once this endpoint is located, however, the software changes the maximum increment to 0.3 ml so that the second endpoint may be quickly located.

The final parameter entered specifies how many consecutive negative second derivatives must be calculated before an endpoint is recognized. This safeguards against the possibility of a false endpoint being detected due to the instability of electrode potential often encountered in the region of the endpoint. As indicated in Figure 2, at the end of the titration a report listing the number of endpoints detected, the endpoint volumes and concentrations calculated is printed on the teletype (TTY).

Upon completion of a titration, the user may have the titration curve plotted on the X-Y recorder. This is accomplished by overlaying the CONVERS plot routine. The plot of the data obtained from the titration shown in Figure 2 is shown in Figure 3. The sample number entered for the titration is reproduced on the plot. In addition, a legend for the plot is printed on the TTY. The numerical coordinates in Figures 3 and 6 were hand lettered.

Titration carried out include titrations of potassium hydrogen phthalate, acetic acid and hydrochloric acids with sodium hydroxide. A Thomas glass electrode (4090-B15) and saturated calomel reference electrode were used.

For ISE calibration experiments, the digital buret is used to create a series of user specified solutions by titrating a stock solution into a blank consisting of either water or a salt solution, whose ionic strength has been adjusted. The software consists of two programs. In the first program (SOLUTION-GEN) the user is asked to input a sample number, the number of electrodes to be calibrated, the blank volume (ml), the highest and lowest concentrations desired (molarity), the change in $-\log$ concentration desired, the stock solution concentration (molarity) and the minimum and maximum volume increments to be added (pulses). The computer uses this information to calculate how many times the buret must be pulsed to create each of the user specified solutions. This data is printed on the TTY and stored in a memory array. A typical computer printout for SOLUTION-GEN is shown in Figure 4.

The second program (CAL-ISE) may then be overlayed. CAL-ISE uses the data stored by SOLUTION-GEN to calibrate the ISEs. The user is asked to input the maximum tolerable electrode potential drift rate (mv/min), the number of minutes to wait for electrode drift to decay below this specified rate, and a delay time (sec) corresponding to the time the computer is to wait between switching an electrode to the meter and sampling. A simplified flow chart of the

calibration process for two electrodes is shown in Figure 5. As can be seen from Figure 5, sampling continues for each electrode until all electrodes produce two consecutive readings that differ by less than the specified drift rate or until the specified maximum time has elapsed.

At the end of the calibration, the $-\log$ concentration-potential data for each electrode is printed on the TTY. In addition, the user is notified during the calibration of those solutions for which satisfactory potential readings could not be obtained. The user may then have the calibration curves plotted on the X-Y recorder by overlaying the plot program.

Calibration data was generated for coated wire ion selective electrodes (CWISEs) of the NO_3^- and ClO_4^- type and for the Orion F^- electrode. The CWISEs were prepared as described elsewhere (11, 12). Potentials were referred to the double junction reference. The outer filling solutions were 0.1M NH_4NO_3 and 0.1M KCl for the ClO_4^- and NO_3^- CWISEs, respectively, and 0.1M KCl for the F^- ISE. Thermostated ($25^\circ\text{C} \pm 0.1^\circ\text{C}$) 300 mL beakers were used as sample cells. The beakers were covered during calibration to prevent evaporative losses.

RESULTS AND DISCUSSION

Meter Response. The 801A uses an integrating analog to digital converter (A/D). This type of A/D has excellent noise rejection, which accounts for the great stability of the digital display, however, ADC is relatively slow. Speed of ADC is not an extremely important consideration here, however, in that ISEs are inherently slow devices. The ADC frequency of the 801A was found to be around 5.5 Hz. Use of the meter print gate as a sampling flag, as described in the experimental section, allows for synchronized data transfer from meter to computer at near this frequency.

The noise present after analog to digital conversion was evaluated by

determination of the standard deviation for both single and averaged readings. The standard deviation for a single reading was found to be around ± 0.03 mv. The standard deviation for a 25 point average was found to be around ± 0.01 mv. Increasing the number of points averaged to 50 did not further lower the standard deviation. The decision to use either single point sampling of meter response or an average of multiple points should be made on the basis of precision desired and electrode characteristics. Changes in the sampling routine may be implemented through minor software modifications. For the calibration curve experiments, a 25 point average was used. For titrations, a single point sample was sufficient.

Titration Experiments

Results of various acid base titrations are summarized in Table I. Each result represents an average of at least four determinations. As indicated, for all titrations relative precisions of under 0.1% were obtained using the automated titration system. Results of titrations of a mixture of acetic and hydrochloric acids are shown in Table II. Good accuracy and precision was obtained for both the HCl and total acid endpoints. This data demonstrates the power of the endpoint locating algorithm when, as shown in Figure 3, the break in a titration curve (HCl endpoint) is not very pronounced.

Calibration Experiments

Calibration curves obtained with the automated system for all electrodes tested were identical to those obtained manually in terms of slope, range of linear response and standard deviation. Calibration curves obtained for four ClO_4^- CWISEs, are shown in Figure 6. In this case, solutions ranging from 10^{-4} to 10^{-2} M ClO_4^- were created. Water was used as the blank solution.

Due to the limitations in resolution and volume of the digital buret, the calibration curves generated by this method have a maximum concentration

dynamic range of about three orders of magnitude. By varying either the standard concentration or the blank volume, however, any concentration region may be investigated. Hence if the response characteristics of an electrode are desired over a larger concentration range, multiple calibrations, with each examining a different portion of the response curve, may be run.

Cost

The approximate costs of the various components used in assembling this system are shown in Table 3. It is important to note that the cost of the computer and associated hardware is less than one-third of the total system cost. Had a minicomputer with comparable hardware and software capabilities been used, the computer hardware cost would have, at least, tripled and the total system cost doubled. It is also important to note that no software expenses have been included (the CONVERS software system may be obtained free of charge).

CONCLUSIONS

A powerful microcomputer-based potentiometric analysis system has been constructed. This system combines the low cost and ease of construction advantages associated with microcomputer-based instrumentation with the high level software capabilities generally only associated with minicomputers. It is the first low cost/high power potentiometric analysis system to be reported and as such should be of great interest to workers involved in both the fundamental and applied aspects of ISE technology.

The ability of the system to locate endpoints in potentiometric titrations and to generate calibration curve data for multiple ISEs has been demonstrated. Other applications such as pK_a and metal-ligand formation constant determination, standard addition and linear (Gran) titrations and automated determination

of ISE selectivity coefficients could easily be implemented. These extensions would merely entail writing the appropriate CONVERS programs. Furthermore, additional instrumentations, such as carousel type sampling devices or equipment associated with flowing stream potentiometry (2), could easily be added to the system. Such hardware extensions would merely entail purchasing additional I/O boards and connecting the desired devices to these boards. Circuit diagrams of interfaces and source listings of programs are available on request.

CONCLUSIONS

A powerful microcomputer-based potentiometric analysis system has been constructed. This system combines the low cost and ease of construction advantages associated with microcomputer-based instrumentation with the high level software capabilities generally only associated with minicomputers. It is the first low cost, high power potentiometric analysis system to be reported and as such should be of great interest to workers involved in both the fundamental and applied aspects of ISE technology.

The ability of the system to locate endpoints in potentiometric titrations and to generate calibration curve data for multiple ISEs has been demonstrated. Other applications such as pH and metal-ligand formation constant determination, standard addition and linear (Gran) titrations and automated determination

REFERENCES

1. D. J. Leggett, Anal. Chem., 50, 718 (1978).
2. P. D. Gaarenstroom, J. C. English, S. P. Perone, and J. W. Bixler, Anal. Chem., 50, 811 (1978).
3. T. F. Christiansen, J. E. Busch, and S. C. Krough, Anal. Chem., 48, 1051 (1976).
4. J. M. Ariano and W. F. Gutknecht, Anal. Chem., 48, 281 (1976).
5. J. W. Frazer, W. Selig, and L. P. Rigdon, Anal. Chem., 49, 1250 (1977).
6. L. P. Rigdon, G. J. Moody, and J. W. Frazer, Anal. Chem., 50, 465 (1978).
7. J. W. Frazer, A. M. Kray, W. Selig, and R. Lim, Anal. Chem., 47, 869 (1975).
8. J. J. Zipper, B. Fleet, and S. P. Perone, Anal. Chem., 46, 2111 (1974).
9. I. Sekerka and J. F. Lechner, Anal. Lett., 7, 463 (1974).
10. Copies of CONVERS software system may be obtained from S. B. Tildon and M. B. Denton, University of Arizona, Tucson, Arizona 85721.
11. H. James, G. Carmack, and H. Freiser, Anal. Chem., 44, 856 (1972).
12. B. M. Kneebone and H. Freiser, Anal. Chem., 45, 449 (1973).

ACKNOWLEDGEMENTS

The authors wish to thank Scot B. Tildon and M. Bonner Denton for valuable consultation on both the software and hardware aspects of this work.

This work was conducted with financial assistance from the Office of Naval Research.




Figure 1. Schematic diagram of instrumentation.

SAMPLE NO. =40
 IF TITRANT CONC KNOWN & SAMPLE FROM SOLUTION X=0
 IF TITRANT CONC UNKNOWN & SAMPLE FROM SOLUTION X=1
 IF TITRANT CONC UNKNOWN & SAMPLE WEIGHED X=2 , X=0
 KNOWN CONC=.10206
 SAMPLE VOL=5
 GO X POINTS PAST ENDPOINT X =15
 MAX DRIFT IN 6 SECS=.01
 QUIT AFTER X TIMES , X=10
 MIN PULSES=10
 MAX PULSES =70
 ENDPT WHEN DDI CHANGES X TIMES , X=4
 TYPE 1 WHEN READY 1
 TOO MUCH DRIFT , READING 00051

 NO. ENDPOINTS=000002
 ENDPOINT VOLS. ARE
 +.193376E+01 +.482602E+01
 UNKNOWN CONCS. ARE
 +.394720E-01 +.985088E-01

Figure 2. Typical computer printout for program TITRATE.




Figure 3. Computer generated plot of titration curve.

Figure 3. Typical computer printout for program TITRATE.




Figure 4. Typical computer printout for program SOLUTION-GEN.




Figure 5. Simplified flowchart for program CAL-ISE.

Figure 6. Computer generated plot of calibration curves for four ClO_4^- coated wire ion selective electrodes.




Table II. Precision and accuracy of automated titrations of acetic and hydrochloric acid mixtures with NaOH.

Relative Error	Relative Precision	
0.10%	0.20%	HCl endpoint
-0.11%	0.13%	Total acid endpoint

Table I. Precision of automated titrations of various acids with 0.1 M NaOH.

Acid (0.1 M)	Relative Precision
Potassium acid phthalate	0.07%
HCl	0.09%
Acetic acid	0.07%

Table II. Precision and accuracy of automated titrations of acetic and hydrochloric acid mixture with NaOH.

	<u>Relative Precision</u>	<u>Relative Error</u>
HCl endpoint	0.20%	0.10%
Total acid endpoint	0.13%	-0.11%

Table III. Approximate costs of computer hardware and peripherals.

<u>ITEM</u>	<u>COST</u>	<u>TOTALS</u>
8080 microcomputer kit	\$ 700.00	
Disk drive	880.00	
Memory boards (3)	420.00	
I/O boards (2)	<u>400.00</u>	
COMPUTER HARDWARE TOTAL		\$2,400.00
Burette	1,500.00	
Meter	1,200.00	
Multiplexer	300.00	
X-Y Recorder	1,500.00	
TTY	<u>1,000.00</u>	
PERIPHERAL TOTAL		\$5,500.00
TOTAL SYSTEM		\$7,900.00

3

OCT 13 1978

472:CAN:716:2
78u472-608

TECHNICAL REPORT DISTRIBUTION LIST, CEN

	<u>No.</u> <u>Copies</u>		<u>No.</u> <u>Copies</u>
Office of Naval Research 800 North Quincy Street Arlington, Virginia 22217 Attn: Code 472	2	Defense Documentation Center Building 5, Cameron Station Alexandria, Virginia 22314	12
ONR Branch Office 536 S. Clark Street Chicago, Illinois 60605 Attn: Dr. George Sandoz	1	U.S. Army Research Office P.O. Box 1211 Research Triangle Park, N.C. 27709 Attn: CRD-AA-IP	1
ONR Branch Office 715 Broadway New York, New York 10003 Attn: Scientific Dept.	1	Naval Ocean Systems Center San Diego, California 92152 Attn: Mr. Joe McCartney	1
ONR Branch Office 1030 East Green Street Pasadena, California 91106 Attn: Dr. R. J. Marcus	1	Naval Weapons Center China Lake, California 93555 Attn: Dr. A. B. Amster Chemistry Division	1
ONR Area Office One Hallidie Plaza, Suite 601 San Francisco, California 94102 Attn: Dr. P. A. Miller	1	Naval Civil Engineering Laboratory Port Hueneme, California 93401 Attn: Dr. R. W. Drisko	1
ONR Branch Office Building 114, Section D 666 Summer Street Boston, Massachusetts 02210 Attn: Dr. L. H. Peebles	1	Professor K. E. Woehler Department of Physics & Chemistry Naval Postgraduate School Monterey, California 93940	1
Director, Naval Research Laboratory Washington, D.C. 20390 Attn: Code 6100	1	Dr. A. L. Slafkosky Scientific Advisor Commandant of the Marine Corps (Code RD-1) Washington, D.C. 20380	1
The Assistant Secretary of the Navy (R,E&S) Department of the Navy Room 4E736, Pentagon Washington, D.C. 20350	1	Office of Naval Research 800 N. Quincy Street Arlington, Virginia 22217 Attn: Dr. Richard S. Miller	
Commander, Naval Air Systems Command Department of the Navy Washington, D.C. 20360 Attn: Code 310C (H. Rosenwasser)	1	Naval Ship Research and Development Center Annapolis, Maryland 21401 Attn: Dr. G. Bosmajian Applied Chemistry Division	
		Naval Ocean Systems Center San Diego, California 91232 Attn: Dr. S. Yamamoto, Marine Sciences Division	

Encl 1

TECHNICAL REPORT DISTRIBUTION LIST, 051C

	<u>No.</u> <u>Copies</u>		<u>No.</u> <u>Copies</u>
Dr. M. B. Denton University of Arizona Department of Chemistry Tucson, Arizona 85721	1	Dr. K. Wilson University of California, San Diego Department of Chemistry La Jolla, California	
Dr. R. A. Osteryoung Colorado State University Department of Chemistry Fort Collins, Colorado 80521	1	Dr. A. Zirino Naval Undersea Center San Diego, California 92132	
Dr. B. R. Kowalski University of Washington Department of Chemistry Seattle, Washington 98105	1	Dr. John Duffin United States Naval Postgraduate School Monterey, California 93940	
Dr. S. P. Perone Purdue University Department of Chemistry Lafayette, Indiana 47907	1	Dr. C. M. Hieftje Department of Chemistry Indiana University Bloomington, Indiana 47401	
		Dr. Victor L. Rehn Naval Weapons Center Code 3813 China Lake, California 93555	
Dr. D. L. Venezky Naval Research Laboratory Code 6130 Washington, D.C. 20375	1	Dr. Christie G. Enke Michigan State University Department of Chemistry East Lansing, Michigan 48824	
Dr. H. Freiser University of Arizona Department of Chemistry Tucson, Arizona 85721		Dr. Kent Eisentraut, MBT Air Force Materials Laboratory Wright-Patterson AFB, Ohio 45433	
Dr. Fred Saalfeld Naval Research Laboratory Code 6110 Washington, D.C. 20375	1	Walter G. Cox, Code 3632 Naval Underwater Systems Center Building 148 Newport, Rhode Island 02840	
Dr. E. Chernoff Massachusetts Institute of Technology Department of Mathematics Cambridge, Massachusetts 02139	1		